

# CLEARIMAGE II® HAND MAGNIFIERS

---

## GENERAL

In response to the need by the low vision population for more cosmetically pleasing and easy to use devices, Designs for Vision, Inc. has designed the new ClearImage II® Hand Magnifiers.

As with the spectacle mounted ClearImage II®, the hand held magnifiers have a field view of 40 mm. The crisp edge to edge clarity enables the patient to read easier and faster. With a choice of two ergonomically correct handles, the patient receives the highest quality reading microscope available.

## FITTING

When ordering, specify power required and the type handle needed:

#1 STRAIGHT or #2 CURVED

## POWERS AVAILABLE

2X - 8X



ClearImage II® Hand Magnifiers are available in your choice of two ergonomically correct handles, CURVED (top) and STRAIGHT (bottom), for individual comfort and ease of reading.



**DESIGNS FOR VISION, INC.**  
1-800-345-4009