

N95 Mask Shape Recommendations

When choosing an N95 mask, use a mask that best accommodates eyewear. The following examples illustrate what to look for in the shape of the mask.



N95 3M 8210




Gerson 2130W



N95 3M 9210/9211



Honeywell ONE-Fit NBW 95V/
Honeywell ONE-Fit N95



Moldex 1500

DVI Sport Nose Pad Adjustments



To make clearance for the mask
 Push **outwards** on the nose pads seen by the **blue arrows**
 and **upwards** shown with the **red arrows**.