

N95 Mask Shape Recommendations

When choosing an N95 mask, use a mask that best accommodates eyewear. The following examples illustrate what to look for in the shape of the mask.



N95 3M 8210



Gerson 2130W



N95 3M 9210/9211



Honeywell ONE-Fit NBW 95V/
Honeywell ONE-Fit N95



Moldex 1500

Nike Skylon Ace Nose Pad Adjustments



Make sure the nose pad (area indicated by the **blue arrows**) rests on the narrow section of the mask that fits over your nose. This should provide the proper clearance between the frame and mask.



DESIGNS FOR VISION, INC.®
Loupes and N95 Masks

Rev. 4.20