

## N95 Mask Shape Recommendations

When choosing an N95 mask, use a mask that best accommodates eyewear. The following examples illustrate what to look for in the shape of the mask.



N95 3M 8210



Gerson 2130W



N95 3M 9210/9211




Honeywell ONE-Fit NBW 95V/  
Honeywell ONE-Fit N95



Moldex 1500

## Skyview Nose Pad Adjustments



To make clearance for the mask  
 Push **outwards** on the nose pads seen by the **blue arrows**  
 and **upwards** shown with the **red arrows**.