

CLEARIMAGE II® MICROSCOPES

GENERAL

This cosmetically pleasing optical system has a 40 mm lens providing the widest possible field of vision that is clear, crisp, and distortion-free from edge to edge. The patient's astigmatic prescription may be incorporated into the system, allowing for ease of reading never before possible. A combination of high index glass and a permanent Magnesium-Fluoride Anti-Reflection coating ensures the blackest, sharpest resolution.

FITTING

The ClearImage II® Microscopes are generally prescribed as a monocular aid due to the close working distance. The balance lens can be clear, frosted or blackout. All ClearImage II® Microscopes are available in your choice of lightweight metal frames as well as the May Fulvue and Half Eye frames. Please refer to our [Microscope Frame Selection Guide](#). In addition, the ClearImage II® Microscope is available in diameters of 40, 35, or 32 mm, in any frame of the doctor's or patient's choice. The vertical measurement of the selected frame will determine the diameter of the microscope.

POWERS AVAILABLE

- 2X (8 diopters)
- 3X (12 diopters)
- 4X (16 diopters)
- 5X (20 diopters)
- 6X (24 diopters)
- 7X (28 diopters)
- 8X (32 diopters)



ClearImage II® Microscope
mounted in May Fulvue frame

