

CLEARIMAGE II® HAND MAGNIFIERS

GENERAL

In response to the need by the low vision population for more cosmetically pleasing and easy to use devices, Designs for Vision, Inc. has designed the new ClearImage II® Hand Magnifiers.

As with the spectacle mounted ClearImage II®, the hand held magnifiers have a field view of 40 mm. The crisp edge to edge clarity enables the patient to read easier and faster. With a choice of two ergonomically correct handles, the patient receives the highest quality reading microscope available.

FITTING

When ordering, specify power required and the type handle needed:

#1 STRAIGHT or #2 CURVED

POWERS AVAILABLE

2X - 8X



ClearImage II® Hand Magnifiers are available in your choice of two ergonomically correct handles, CURVED (top) and STRAIGHT (bottom), for individual comfort and ease of reading.



DESIGNS FOR VISION, INC.®
1-800-345-4009