



# CLEARIMAGE II® FRAME SELECTION GUIDE

May Fulvue II Frame  
with 5X ClearImage II®  
Microscope



**NEW** Nadler Mauve Frame  
with 3X ClearImage II®  
Microscope



Gold Wire Half Eye Frame  
with 5X ClearImage II®  
Microscope



Wire 047 Gold Frame  
with 6X ClearImage II®  
Microscope








Yeoman Frame  
with 3X ClearImage II®  
Microscope



May Fulvue II Frame  
with 8X ClearImage II® Telephoto  
Microscope and Focus Tube



# CLEARIMAGE II® FRAME SELECTION GUIDE Sizes and Specifications

<b>May Fulvue II</b> 	<b>COLOR:</b> Tortoise		Adjustable nose pads available upon request				
	<b>FRAME SIZE:</b> 46-22 or 46-24		<b>TEMPLE SIZES:</b> 135, 140, 145, or 150				
	<b>LENS DIMENSIONS:</b> A - 45 mm B - 45 mm ED - 45 mm						
	Included in price of ClearImage II® / \$20.00 when purchased separately						
<b>Gold Wire Half Eye</b> 	<b>COLOR:</b> Gold		<b>NOT available for telephoto microscopes</b>				
	<b>FRAME SIZE:</b> 45-19		<b>TEMPLE SIZES:</b> 145 or 150				
	<b>LENS DIMENSIONS:</b> A - 43 mm B - 27 mm ED - 45 mm						
	Included in price of ClearImage II® / \$35.00 when purchased separately						
<b>Nadler</b> 	<b>COLORS:</b> Copper (top left) Silver (bottom left) Mauve (on front page)						
	<b>FRAME SIZE:</b> 48-17						
	<b>LENS DIMENSIONS:</b> A - 48 mm B - 39 mm ED - 46 mm						
	<b>FRAME SIZE:</b> 51-17						
	<b>LENS DIMENSIONS:</b> A - 50 mm B - 41 mm ED - 49 mm						
	<b>TEMPLE SIZES:</b> 135 or 140						
Additional \$35.00 with ClearImage II® / \$55.00 when purchased separately							
<b>Wire 047</b> 	<b>COLORS:</b> Tortoise (on left) / Gold (on front page)						
	<b>FRAME SIZE:</b> 51-18		<b>TEMPLE SIZE:</b> 140				
	<b>LENS DIMENSIONS:</b> A - 50 mm B - 44 mm ED - 51 mm						
	Additional \$25.00 with ClearImage II® / \$45.00 when purchased separately						
<b>Yeoman</b> 	<b>COLOR:</b> Black						
	<b>AVAILABLE FRAME SIZES:</b>						
	42-20 42-22	44-20 44-22	46-20 46-22 46-24	48-20 48-22 48-24	50-22 50-24	52-22 52-24	54-22 54-24
	<b>TEMPLE SIZES:</b> 5.50 / 5.75 / 6.00 / 6.25 / 6.50 / 6.75 / 7.00						
	The difference between the A (horizontal) and B (vertical) is 6 mm.						
	Additional \$60.00 with ClearImage II® / \$80.00 when purchased separately						

