

SUN ACCESSORIES

Sunglasses are often necessary for the Low Vision patient who is bothered by glare when outside. Although the Low Vision aids themselves can be ordered with a tint, as well as special transition lenses, this is sometimes not enough, especially when the aid is for general use, i.e., both indoors and outside. Designs for Vision, Inc has two different styles of sunglasses that can be used with either telescopes or microscopes. These are the Slip-Behind and the Wrap-Around. Besides providing glare reduction, the Wrap-Around has the additional advantage of covering the aid – a definite benefit for those patients who are very cosmetic conscious.

Slip-Behind Sunglasses: Designed to go behind the frame, these sunglasses also have the feature of a wrap-around edge which helps to reduce the glare coming in from the side of the patient.

Wrap-Around Sunglasses: These are appropriate for any prescription. When used with a telescopic aid, a hole can be drilled in the lens, through which the telescopic unit is placed (must be specified on order if desired). This allows the sunglass to fit flush against the frame that the patient is wearing. This sunglass also has a side piece which helps to reduce glare coming in from the side of the patient.

COLORS AVAILABLE

Slip-Behind Sunglasses - grey/green (standard), amber, yellow

Wrap-Around Sunglasses - grey/green (standard), amber

Other colors may be special ordered for both.



Slip-Behind Sunglasses



Wrap-Around Sunglasses
Custom Drilled
to Fit over Telescopes



DESIGNS FOR VISION, INC.
1-800-345-4009