



TELESCOPE FRAME SELECTION GUIDE

NEW Nadler Copper Frame
with Chocolate 2.2X
Bioptic I Telescope




Rimless Silver Frame
with Silver 2.2X
Bioptic I Telescope



Titanium 077 Chocolate
Chrome Frame
with Chocolate 3.0X
Bioptic Telescope



Titanium 109DB Chocolate
Chrome Frame
with Chocolate 2.2X
Bioptic I Telescope



Yeoman Black Frame
with Black 2.2X
Bioptic I Telescope



Wire 047 Gold Frame
with Black 3.0X
Bioptic Telescope










Titanium 077
Silver Chrome Frame
with Beige 2.2X
Full Diameter Telescope



Titanium 062
Chocolate Chrome Frame
with Chocolate 2.2X
Bioptic I Telescope



TELESCOPE FRAME SELECTION GUIDE Sizes and Specifications

Yeoman		COLOR: Black						
		AVAILABLE FRAME SIZES:						
		42-20 42-22	44-20 44-22	46-20 46-22 46-24	48-20 48-22 48-24	50-22 50-24	52-22 52-24	54-22 54-24
		TEMPLE SIZES: 5.50 / 5.75 / 6.00 / 6.25 / 6.50 / 6.75 / 7.00						
		The difference between the A (horizontal) and B (vertical) is 6 mm.						
Included in price of Telescopic System / \$80.00 when purchased separately								
Nadler		COLORS: Silver (top left) / Mauve (bottom left) / Copper (on front page)						
		FRAME SIZE: 48-17						
		LENS DIMENSIONS: A - 48 mm B - 39 mm ED - 46 mm						
		FRAME SIZE: 51-17						
		LENS DIMENSIONS: A - 50 mm B - 41 mm ED - 49 mm						
TEMPLE SIZES: 135 or 140								
Additional \$25.00 with Telescopic System / \$55.00 when purchased separately								
Wire 047		COLORS: Tortoise (on left) / Gold (on front page)						
		FRAME SIZE: 51-18	TEMPLE SIZE: 140					
		LENS DIMENSIONS: A - 50 mm B - 44 mm ED - 51 mm						
		Additional \$15.00 with Telescopic System / \$45.00 when purchased separately						
Titanium 062		COLOR: Chocolate Chrome						
		FRAME SIZE: 54-19	TEMPLE SIZES: 140 or 150					
		LENS DIMENSIONS: A - 54 mm B - 44 mm ED - 56 mm						
		Additional \$70.00 with Telescopic System / \$100.00 when purchased separately						
Wire 109DB:		COLORS: Wire Frames: Silver / Blue / Rose / Gold (on left) Titanium Frame: Chocolate (on front pg)						
		NOT available for Expanded Field Telescopes						
		FRAME SIZE: 52-17	TEMPLE SIZES: 135 or 140					
		LENS DIMENSIONS: A - 52 mm B - 35 mm ED - 53 mm						
		Wire Frames: Additional \$20.00 with Telescopic System \$50.00 when purchased separately Titanium Frames: Additional \$70.00 with Telescopic System \$100.00 when purchased separately						
Titanium 077		COLOR: Chocolate Chrome (left) Silver Chrome (on front page)						
		NOT available for Expanded Field Telescopes						
		FRAME SIZE: 48-18	TEMPLE SIZE: 135 or 140					
		LENS DIMENSIONS: A - 48 mm B - 33 mm ED - 47 mm						
Additional \$70.00 with Telescopic System / \$100.00 when purchased separately								
Rimless		COLORS: Antique Gold (top left) / Black (bottom left) / Silver (on front page)						
		TEMPLE SIZES: 135, 140, 145 or 150						
		LENS DIMENSIONS LoRusso: A-52mm B-36mm ED-56mm / Shelby: A-51mm B-38mm ED-51mm						
		Plano additional \$110.00 or Single Vision additional \$135 with Telescopic System \$200.00 when purchased separately and includes A/R coated Trivex Lenses						

