



# TELESCOPE FRAME SELECTION GUIDE

**NEW Nadler Copper Frame**  
with Chocolate 2.2X  
Bioptic I Telescope



**NEW Steam Purple Frame**  
with **NEW** Black Bio Micro EF  
4.5X expanded field telescope



**NEW Steam Black Frame**  
with Black 3.0X  
Bioptic Telescope



**NEW Botta Chocolate  
Titanium Frame**  
with Chocolate 2.2X  
Bioptic I Telescope



**Yeoman Black Frame**  
with Black 2.2X  
Bioptic I Telescope



**Wire 047 Gold Frame**  
with Black 3.0X  
Bioptic Telescope










**Titanium 077  
Silver Chrome Frame**  
with Beige 2.2X  
Full Diameter Telescope



**Titanium 062  
Chocolate Chrome Frame**  
with Chocolate 2.2X  
Bioptic I Telescope



# TELESCOPE FRAME SELECTION GUIDE Sizes and Specifications

|   |   |  |   |                         |                         |                |                |                |
|---|---|--|---|-------------------------|-------------------------|----------------|----------------|----------------|
| Yeoman  |    | <b>COLOR:</b> Black  |   |                         |                         |                |                |                |
|   |   | <b>AVAILABLE FRAME SIZES:</b>  |   |                         |                         |                |                |                |
|   |   | 42-20<br>42-22   | 44-20<br>44-22  | 46-20<br>46-22<br>46-24 | 48-20<br>48-22<br>48-24 | 50-22<br>50-24 | 52-22<br>52-24 | 54-22<br>54-24 |
|   |   | <b>TEMPLE SIZES:</b> 5.50 / 5.75 / 6.00 / 6.25 / 6.50 / 6.75 / 7.00                  |   |                         |                         |                |                |                |
| The difference between the A (horizontal) and B (vertical) is 6 mm. |   |  |   |                         |                         |                |                |                |
| Nadler  |    | <b>COLORS:</b> Silver, Mauve, Copper (on front page)                                 |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 48-17 <b>LENS DIMENSIONS:</b> A - 48 mm B - 39 mm ED - 46 mm      |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 51-17 <b>LENS DIMENSIONS:</b> A - 50 mm B - 41 mm ED - 49 mm      |   |                         |                         |                |                |                |
|   |   | <b>TEMPLE SIZES:</b> 135 or 140  |   |                         |                         |                |                |                |
| Steam   |    | <b>COLORS:</b> Brown, Gray, Navy, Black (on front page) and Purple (on front page)   |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 52-16   |   |                         |                         |                |                |                |
|   |   | <b>TEMPLE SIZE:</b> 137  |   |                         |                         |                |                |                |
|   |   | <b>LENS DIMENSIONS:</b> A - 52 mm B - 42 mm ED - 52 mm                               |   |                         |                         |                |                |                |
| Wire 047  |    | <b>COLORS:</b> Tortoise (on left) / Gold (on front page)                             |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 51-18   | <b>TEMPLE SIZE:</b> Gold 140 only, Tortoise: 135, 140, 145 or 150 |                         |                         |                |                |                |
|   |   | <b>LENS DIMENSIONS:</b> A - 50 mm B - 44 mm ED - 51 mm                               |   |                         |                         |                |                |                |
| Titanium 062  |   | <b>COLOR:</b> Chocolate Chrome   |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 54-19   | <b>TEMPLE SIZES:</b> 140 or 150                                   |                         |                         |                |                |                |
|   |   | <b>LENS DIMENSIONS:</b> A - 54 mm B - 44 mm ED - 56 mm                               |   |                         |                         |                |                |                |
| Botta   |   | <b>COLORS:</b><br>Wire Frames: Silver, Blue, Gold<br>Titanium Frame: Chocolate, Rose |   |                         |                         |                |                |                |
|   |   | <b>NOT available for Expanded Field Telescopes</b>                                   |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 52-17   | <b>TEMPLE SIZES:</b> 135 or 140                                   |                         |                         |                |                |                |
|   |   | <b>LENS DIMENSIONS:</b> A - 52 mm B - 35 mm ED - 55 mm                               |   |                         |                         |                |                |                |
| Titanium 077  |  | <b>COLOR:</b> Chocolate Chrome (left)<br>Silver Chrome (on front page)               |   |                         |                         |                |                |                |
|   |   | <b>FRAME SIZE:</b> 48-18   | <b>TEMPLE SIZE:</b> 135 or 140                                    |                         |                         |                |                |                |
|   |   | <b>LENS DIMENSIONS:</b> A - 48 mm B - 33 mm ED - 47 mm                               |   |                         |                         |                |                |                |



**Designs for Vision, Inc.**

4000 Veterans Memorial Hwy. • Bohemia, NY 11716-1024 • 1.631.585.3300

1.800.345.4009 • FAX 1.631.585.5704 • [LowVision@dvimail.com](mailto:LowVision@dvimail.com)

[DesignsForVision.com](http://DesignsForVision.com)