



TELESCOPE FRAME SELECTION GUIDE

NEW Nadler Copper Frame
with Chocolate 2.2X
Bioptic I Telescope



NEW Steam Purple Frame
with **NEW** Black Bio Micro EF
4.5X expanded field telescope



NEW Steam Black Frame
with Black 3.0X
Bioptic Telescope



Botta Chocolate
Titanium Frame
with Chocolate 2.2X
Bioptic I Telescope



Yeoman Black Frame
with Black 2.2X
Bioptic I Telescope



Wire 047 Gold Frame
with Black 3.0X
Bioptic Telescope










Titanium 077
Silver Chrome Frame
with Beige 2.2X
Full Diameter Telescope



Titanium 062
Chocolate Chrome Frame
with Chocolate 2.2X
Bioptic I Telescope



TELESCOPE FRAME SELECTION GUIDE Sizes and Specifications

Yeoman		COLOR: Black						
		AVAILABLE FRAME SIZES:						
		42-20 42-22	44-20 44-22	46-20 46-22 46-24	48-20 48-22 48-24	50-22 50-24	52-22 52-24	54-22 54-24
		TEMPLE SIZES: 5.50 / 5.75 / 6.00 / 6.25 / 6.50 / 6.75 / 7.00						
		The difference between the A (horizontal) and B (vertical) is 6 mm.						
Included in price of Telescopic System / <i>\$80.00 when purchased separately</i>								
Nadler		COLORS: Silver, Mauve, Copper (on front page)						
		FRAME SIZE: 48-17 LENS DIMENSIONS: A - 48 mm B - 39 mm ED - 46 mm						
		FRAME SIZE: 51-17 LENS DIMENSIONS: A - 50 mm B - 41 mm ED - 49 mm						
		TEMPLE SIZES: 135 or 140						
		Additional \$25.00 with Telescopic System / <i>\$55.00 when purchased separately</i>						
Steam		COLORS: Brown, Gray, Navy, Black (on front page) and Purple (on front page)						
		FRAME SIZE: 52-16						
		TEMPLE SIZE: 137						
		LENS DIMENSIONS: A - 52 mm B - 42 mm ED - 52 mm						
		Additional \$70.00 with Telescopic System / <i>\$100.00 when purchased separately</i>						
Wire 047		COLORS: Tortoise (on left) / Gold (on front page)						
		FRAME SIZE: 51-18	TEMPLE SIZE: Gold 140 only, Tortoise: 135, 140, 145 or 150					
		LENS DIMENSIONS: A - 50 mm B - 44 mm ED - 51 mm						
		Additional \$15.00 with Telescopic System / <i>\$45.00 when purchased separately</i>						
Titanium 062		COLOR: Chocolate Chrome						
		FRAME SIZE: 54-19	TEMPLE SIZES: 140 or 150					
		LENS DIMENSIONS: A - 54 mm B - 44 mm ED - 56 mm						
		Additional \$70.00 with Telescopic System / <i>\$100.00 when purchased separately</i>						
Botta		COLORS:						
		Wire Frames: Silver, Blue, Gold Titanium Frame: Chocolate, Rose	NOT available for Expanded Field Telescopes					
		FRAME SIZE: 52-17	TEMPLE SIZES: 135 or 140					
		LENS DIMENSIONS: A - 52 mm B - 35 mm ED - 55 mm						
		Wire Frames: Additional \$20.00 with Telescopic System <i>\$50.00 when purchased separately</i> Titanium Frames: Additional \$70.00 with Telescopic System <i>\$100.00 when purchased separately</i>						
Titanium 077		COLOR: Chocolate Chrome (left) Silver Chrome (on front page)						
		FRAME SIZE: 48-18	TEMPLE SIZE: 135 or 140					
		LENS DIMENSIONS: A - 48 mm B - 33 mm ED - 47 mm						
		Additional \$70.00 with Telescopic System / <i>\$100.00 when purchased separately</i>						



Designs for Vision, Inc.

4000 Veterans Memorial Hwy. • Bohemia, NY 11716-1024 • 1.631.585.3300

1.800.345.4009 • FAX 1.631.585.5704 • LowVision@dvimail.com

DesignsForVision.com